



Dealer / Location:		T:	. F:
DO# / Clients	Contact	Emaile	

TRAKZ TUB SLIDERS CONSIST OF THREE MAJOR COMPONENTS:

Mobile Bases – Include height-adjustable aluminum frame, soft contoured seat cushion, molded vinyl back upholstery, swing-away removable padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. The Trakz-AT Mobile Base also includes a rotational tilt system with a range of 0° to 35° and a head support.

Tub Base – Includes a toolless, height-adjustable aluminum frame with horizontal stabilizers.

Bridge – The patented Zbrilok safety system prevents unintentional seat release from either base.

Zbrilok™

The patented latching safety system, Zbrilok, ensures that the sliding seat assembly cannot be moved onto the bridge, from either the tub base or mobile base, unless the bridge is fully locked to both bases. This prevents the seat assembly from being released inadvertantly and potentially causing injury to the client. The Raz Trakz is the only tub slider that has this safety feature.


1. Models (NOTE: must select either T188 AND T1BASE or T388 AND T3BASE)

T188	Trakz-AP (Attendant Propel) Non-Tilting Tub Slider Mobile Base	\$2,895
T1BASE	Tub Base for Trakz-AP includes adjustable-height frame and horizontal stabilizers (without tools)	\$2,595
T1BASEV	Tub Base for Trakz-AP with Vertical Stabilizers (Tub Base for Step-in-Shower)	\$2,595
T388	Trakz-AT (Attendant Tilt) Tilting Tub Slider Mobile Base with 35° of rotational tilt and head support	\$4,125
T3BASE	Tub Base for Trakz-AT includes adjustable-height frame and horizontal stabilizers (without tools)	\$2,595
T3BASEV	Tub Base for Trakz-AT with Vertical Stabilizers (Tub Base for Step-in-Shower)	\$2,595

2. Slide Orientation: Must select one

TSOL	Slide Orientation - Left	No charge
TSOR	Slide Orientation - Right	No charge

Changing a Trakz-AT from right to left or left to right is easy. The tilt triggers and tilt assist handle are the only components that need to be relocated.



Slide Orientation – Left



Date:____

Slide Orientation – Right

3. Seats


Visco-Foam Interface Seats (Increases the seat-to-floor height by 1.75") Only available for Trakz-AT				
TVFI17FRO	18"W x 17"D – Front Access	\$475		
TVFI17BRI	18"W x 17"D - Front Bridged	\$475		
TVFI17RRO	18"W x 17"D – Rear Access	\$475		
TVFI20FRO	18"W x 20"D – Front Access	\$545		
TVFI20BRI	18"W x 20"D – Front Bridged	\$545		
TVFI20RRO	18"W x 20"D – Rear Access	\$545		
Trakz Custom Vi	Trakz Custom Visco Foam Seat - Similar to Cat Seat with oblong aperture available 18" x 17" or 18" x 18"			
TZCSEAT		\$1145		
NOTE: ** will incre	NOTE: ** will increase STF by 1"			

RAZ TRAZK - 23-01 -25-CP

3. Bridge Lengths (See next page for bridge length calculation)

TBL07	Bridge length – 7"			\$490
TBL08	Bridge length – 8"			\$490
TBL09	Bridge length – 9"			\$490
TBL10	Bridge length – 10"			\$490
TBL12	Bridge length – 12"			\$490
TBL14	Bridge length – 14"			\$490
TBL16	Bridge length – 16"			\$490
TBLC	Bridge length – Custom	Enter length required"		\$825

Bridge Length without toilet interference



Calculated Bridge Length (CBL) – without toilet bowl interference

Minimum CBL = $D_1 + 2" - 10.25"$


Example: If $D_1 = 15$ "

Minimum CBL = 15 + 2" - 10.25" = 6.75" Select: TBL07 Bridge length – 7" or longer D₁ = Distance between center of tub and outside of tub wall

10.25" = Half the width of the tub base

BL = Bridge Length (round up to the nearest available size [See Sec 2 - previous page]

Bridge Length with toilet interference



Calculated Bridge Length (CBL) – with toilet bowl interference

Minimum CBL = $D_2 - 10.25'' - 3.5''$

Example: If $D_2 = 25''$ CBL = 25'' - 10.25'' - 3.5''

CBL = 11.25''

Select: TBL11 Bridge length -11.5"

BL = Bridge Length (round up to the nearest available size) [See Sec 2 - previous page]

D₂ = Distance between center of tub and near side of toilet bowl

4. Trakz Back Support Options

The standard Tra	kz back support includes 18"W molded, vinyl upholstery, wrap-around and removable padded	d arm supports.
TFB1820	Fixed Back with molded vinyl upholstery – 18"W x 20"H	Standard
TFBU1820	Fixed Back with Tension-Adjustable Upholstery – 18"W x 20"H	\$195
TFBSB1820	Fixed Back with Symphony Back – 18"W x 20"H	\$315

5. Arm Supports – Wrap-around, swing-away, removable

TZASWL	Arm Supports (pr) Wraparound – Long	Standard
TZASWS	Arm Supports (pr) Wraparound – Short	N/C
Arm Support	Pads	
TZARMLL	Molded pad long 13" – Left (fits TZASWL)	Standard
TZARMLR	Molded pad long 13" – Right (fits TZASWL)	Standard
TZARMSL	Molded pad short ¹ 10" – Left (fits TZASWS)	N/C
TZARMSR	Molded pad short ¹ 10" – Right (fits TZASWS)	N/C
TZTAPLL	Tubular pad long 13" – Left (in lieu of molded pads, fits TZASWL)	N/C
TZTAPLR	Tubular pad long 13" – Right(in lieu of molded pads, fits TZASWL)	N/C
TZTAPSL	Tubular pad short 10" – Left (in lieu of molded pads, fits TZASWS)	N/C
TZTAPSR	Tubular pad short 10" – Right (in lieu of molded pads, fits TZASWS)	N/C
TZFAPLL	Molded Flat Armpad Long 14"L x 3.25"W – Left (fits TZASWL)	\$91
TZFAPLR	Molded Flat Armpad Long 14"L x 3.25"W – Right (fits TZASWL)	\$91
TZFAPSL	Molded Flat Armpad Short 10"L x 3"W – Left (fits TARMS)	\$91
TZFAPSR	Molded Flat Armpad Short 10"L x 3"W – Right (fits TARMS)	\$91
TZTRL	Molded Arm Trough 18.5"L x 4.5"W –Left (fits TZASWL)	\$125

TZTRR	Molded Arm Trough 18.5"L x 4.5"W – Right (fits TZASWL)	\$125	
TZPAML	Pivoting Arm Mount ² – Left (for TZFAPLL / TZTRL)	\$185	
TZPAMR	Pivoting Arm Mount ² – Right (for TZFAPLR / TZTRR)	\$185	
TZAPSB18	Anterior Posture Support Bar (for ZARML only)	\$245	
	¹ Not compatible with TZFAP / TZTR / TZPAM ² Recommended when Lateral Thoracic Supports (TLTL/S) and Molded Arm Troughs (TZTR) are both ordered		

6. Foot / Leg Supports and Accessories

_		• •	
	TZFSL	Foot Support – Left 16" – 20.5"	Standard
	TZFSR	Foot Support – Right 16" – 20.5"	Standard
	TZFET2L	Footrest Extension Tubes ¹ 2" longer – Left (for TZFSL)	\$31
	TZFET2R	Footrest Extension Tubes ¹ 2" longer – Right (for TZFSR)	\$31
	TZMFX8L	MFX8 Foot Support – Left 8" – 12"	\$57
	TZMFX8R	MFX8 Foot Support – Right 8" – 12"	\$57
	TZMFX12L	MFX12 Foot Supports – Left 12" – 16.5"	\$57
	TZMFX12R	MFX12 Foot Supports – Right 12" – 16.5"	\$57
	TZADML	Adj Angle/Depth Molded Footplates ² – Left (for TZFSL / TZMFX12L)	\$110
	TZADMR	Adj Angle/Depth Molded Footplates ² – Right (for TZFSR / TZMFX12R)	\$110
	TZFET2AL	Footrest Extension Tubes 2" longer – Left (for TZADML)	\$31
	TZFET2AR	Footrest Extension Tubes 2" longer – Right (for TZADMR)	\$31
	TZELRL	Elevating Leg Support – Left (includes TZADML)	\$550
	TZELRR	Elevating Leg Support – Right (includes TZADMR)	\$550
	TZRLSL	Universal Residual Limb Support – Left	\$215
	TZRLSR	Universal Residual Limb Support – Right	\$215
	TZCS18	Calf Strap, Snap-On – Bodypoint Aeromesh	\$150
	TZCSIC18	Calf Strap, Snap-On – Infection Control	\$140
	TZHS1820	H-Strap – hook and loop attachment	\$102
	TZAHL	Adjustable Heel Loops (pr)	\$70
	TZNC	Neoprene Cover for Standard Footplate (Each Per Side)	\$51
	TZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates (Each Per Side)	\$51
1	Not compatible v	with TZMFX8 ² Not compatible with TZMFX8 / TZVS	
_			

7. Positioning Pads

Head Support	S		
THSP1	Head Support Pad for T388 – Small (10" x 4")	Standard	
THSP2	Head Support Pad for T388 – Large (12" x 6")	\$55	
THSA1	Head Support Assembly for T188 – with Small Pad (10" x 4")	\$535	
THSA2	Head Support Assembly for T188 – with Large Pad (12" x 6")	\$590	
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, select ZLAEP extension plate.			
TZLTSWA1L	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Left	\$365	
TZLTSWA1R	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Right	\$365	
TZLTSWA2L	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Left	\$365	
TZLTSWA2R	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Right	\$365	

	ates for Lateral Thoracic Supports allow for an additional 1.5" medially and can be used with Lateral Thoracic Supports	F
TZLAEPL	Lateral Extension Plate – Left	\$27
TZLAEPR	Lateral Extension Plate – Right	\$27

8. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
TZPBM	Pelvic Belt – Medium (2-piece w/ side-release buckle)	\$60	TZCBM	Chest Belt ² – Medium (1-piece w/ hook and loop closure)	\$43
TZPBL	Pelvic Belt – Large (2-piece w/ side-release buckle)	\$55	TZCBL	Chest Belt ² – Large (1-piece w/ hook and loop closure)	\$55
TZPBICM	Infection Control Pelvic Belt ¹ – Medium (2-piece w/ side-release buckle)	\$105	TZCBICM	Infection Control Chest Belt ^{1,2} – Medium (1-piece w/ side-release buckle)	\$85
TZPBICL	Infection Control Pelvic Belt ¹ – Large (2-piece w/ side-release buckle)	\$105	TZCBICL	Infection Control Chest Belt ^{1,2} – Large (1-piece w/ side-release buckle)	\$95
TZPBBPM	Pelvic Belt – Bodypoint Aeromesh – Med	\$285		•	
TZPBBPL	Pelvic Belt – Bodypoint Aeromesh – Large	\$285	TZCBBPM	Chest Belt – Bodypoint Aeromesh – Med	\$285
.2.0012	. c.n.e ben bodypolite/terofficsit Edige	Ţ <u></u>	TZCBBPL	Chest Belt – Bodypoint Aeromesh – Large	\$295

¹ Non-absorbent belts – recommended for infection-control environments

9. Other

TZIHW	Whizard – winged urine deflector – Has patented lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$75
TZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$165
TZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$375

10. Credits and Omissions

TZ905	Omit Commode Pan and Mounting Brackets	-\$20	TZ911	Omit Molded Head Support Pad (T388 only)	-\$20
TZ907	Omit Foot Supports – if no other Foot / Leg Support Option is selected	-\$45	T912	Omit Complete Head Support Assembly (T388 only)	-\$50

Some assembly of Raz products is required.


11. Notes:

² Optional tension – Adjustable Upholstery includes loops to retain a chest belt (ZCBM / ZCBL). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.


Specifications	Trakz-AP Attendant Propel	Trakz-AT Attendant Tilt
Weight Capacity	35	50lbs
Overall Width - with arm supports - without arm supports		23"
Overall Length with standard footrests	35.5"	36"
Tilt Range	N/A	35°
Back Support Height	:	20"
Bridge Length	9.5", 12.5", 16.5	' + custom lengths
Underside of Bridge to Top of Seat	6.25"	8.5"
Underside of Bridge to Floor - high base configuration - low base configuration	· ·	54cm) increments .54cm) increments
Seat Width		18"
Seat Depth	17"	(43cm)
Clearance Height Under Seat - high base configuration - low base configuration * without brackets +0.75"		in 1" increments* in 1" increments*
Arm Support Height Above Seat		10"
Seat-to-Floor Height - high base configuration - low base configuration	22" – 25" in 1" increments 17" – 21" in 1" increments	24.5" – 27.5" in 1" increments 19.5" – 23.5" in 1" increments
Product Weight	77.5 lbs	95 lbs
Frame Materials - mobile base frame - tub base frame - bridge - back frame and footrests - hardware	Alur Alur Stainl	minum minum minum ess Steel ess Steel

RAZ TRAZK - 23-01 -25-CP

Client Information









lmage 2

Left	Right	Image 1	Image 2
		A Buttock / Thigh depth	E Chest width
		B Lower leg length	F Hip width
		C Seat to elbow	
		D Seat to Inferior Scapula	

Common Bathroom Layouts RAZTRAZK-23-01-25-CP Fixtures on common wall




		"
Α	Overall length of bathtub	
В	Length of bathtub floor	"
С	Width of bathtub floor	"
D	Overall width of bathtub	"
Е	Distance between bathtub and toilet bowl edge	"
F	Widest part of toilet bowl	"
G	Distance between edge of toilet and closest wall	"
Н	Distance between edge of toilet and opposite wall	"
I	Distance from center of toilet bowl and center of bathtub	"
J	Distance from top edge of bathtub to bathtub floor	"
K	Distance from top edge of bathtub to bathroom floor	"
L	Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)	"
М	Width of door opening	"




Α	Overall length of bathtub	"
В	Length of bathtub floor	"
С	Width of bathtub floor	"
D	Overall width of bathtub	"
Е	Distance between bathtub and toilet bowl edge	"
F	Widest part of toilet bowl	"
G	Distance between edge of toilet and closest wall	"
Н	Distance between edge of toilet and opposite wall	"
I	Distance from center of toilet bowl and center of bathtub	"
J	Distance from top edge of bathtub to bathtub floor	"
K	Distance from top edge of bathtub to bathroom floor	"
L	Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)	"
М	Width of door opening	"

RAZ TRAZK - 23-01 -25-CP

Common Bathroom Layouts Fixtures on two walls



Α	Overall length of bathtub	"
В	Length of bathtub floor	"
С	Width of bathtub floor	"
D	Overall width of bathtub	"
Е	Distance between bathtub and toilet bowl edge	"
F	Widest part of toilet bowl	"
G	Distance between edge of toilet and closest wall	"
Н	Distance between edge of toilet and vanity	"
ı	Distance from top edge of bathtub to bathtub floor	"
J	Distance from top edge of bathtub to bathroom floor	"
K	Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)	"
L	Width of door opening	"



Α	Overall length of bathtub	ıı ıı
В	Length of bathtub floor	"
С	Width of bathtub floor	"
D	Overall width of bathtub	"
Е	Distance between bathtub and toilet bowl edge	"
F	Widest part of toilet bowl	"
G	Distance between edge of toilet and closest wall	"
Н	Distance between edge of toilet and vanity	"
ı	Distance from top edge of bathtub to bathtub floor	"
J	Distance from top edge of bathtub to bathroom floor	"
К	Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)	"
L	Width of door opening	II.

	RAZ TRAZK - 23-01 -2		
	oom Sket		

Overall length of bathtub	"
Length of bathtub floor	"
Width of bathtub floor	"
Overall width of bathtub	"
Distance between bathtub and toilet bowl edge	"

Wi	dest part of toilet bowl	"
	stance between edge of let and closest wall	"
1	stance between edge of let and opposite wall	"
	stance from center of toilet wl and center of bathtub	"

Distance from top edge of bathtub to bathtub floor	"
Distance from top edge of bathtub to bathroom floor	"
Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)	II.
Width of door opening	"