



Dealer / Location:		T:	. F:
PO# / Client:	Contact:	Email:	

Züm-AP Folding Attendant Propel Mobile Shower Commode Chair. Includes aluminum frame, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. User weight capacity 350 lbs. See our Evaluation and Medical Justification Form.

1. Chair Code

Z11118	Züm-AP w/ 18"W frame – will fit over most toilets	\$2,875
--------	---	---------

2. Züm-AP Travel Kits – Must select 15"H Back Support (ZBZ1815)

ZZUMAPTK	Züm-AP Travel Kit – Includes suitcase, component bags, removable casters. Commode Pan	\$525
	will not fit in Travel Kit suitcase.	

2. Seats - Standard

Molded Seats (t	hickness ~2.5"-3")		Visco-Foam Inter	rface Seats (thickness ~3.5")	
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$475
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$475
ZMS20FRO	18"W x 20"D – Front Access	\$150	ZVFI17LTO	18"W x 17"D – Left Access	\$475
ZMS20BRI	18"W x 20"D – Front Bridged	\$150	ZVFI17RTO	18"W x 17"D – Right Access	\$475
Specialty Seats			ZVFI17RRO	18"W x 17"D – Rear Access	\$475
ZSOM	Shower-Only Seat – Molded	\$440	ZVFI20BRI	18"W x 20"D – Front Bridged	\$545
	18"W x 16"-21"D		ZVFI20FRO	18"W x 20"D – Front Access	\$545
ZEJRM	E&J Replica Seat – Molded	\$215	ZVFI20LTO	18"W x 20"D – Left Access	\$545
– 18"W x 17"D		ZVFI20RTO	18"W x 20"D – Right Access	\$545	
			ZVFI20RRO	18"W x 20"D – Rear Access	\$545

3. Seats - Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

ZCSEAT	Custom Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1145
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are ≤ 21"	\$1305

Custom Specifications (for above custom seats)

Custom Seat Width: _ _ 18" / 20" (cannot be more than 2" wider than frame width)

Custom Seat Depth: __ 16" / 17" / 18" / 19" / 20"

Access Opening: __ FR - Front / LT - Left / RT - Right / RR - Rear / BR - Bridged Access Opening Width: __ 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)

Aperture Width: __ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than 18"W.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support. Consider using Arm Support Spacer Kit (ZASK05/ZASK10) to increase width between the arm supports by 1" or 2".
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.



4. Fixed Back Frame – These are called "Fixed" Backs because the arms supports are not adjustable in height. Inludes molded vinyl back support. Back heights are nominally measured from top of standard molded seat to top of back upholstery.

ZFBZ1817	Fixed Back – 17" Height ¹	Included
¹ Arm Support Lock	Kit not included [See Sec.5]	

5. Arm Support Accessories - Must order Tubular Arm Support Pads or Long / Short Arm Support Pads for compatibility with Zum Travel Kits (ZZUMAPTK)

ZARMS	Flip-Up Arm Supports – 10"L pad (pr) (not compatible with ZFL / ZFS /	Std	ZTRL	Molded Arm Trough – Left (18.5″L x 4.5″W – fits ZARML)	\$125
	ZTR / ZPAM)		ZTRR	Molded Arm Trough – Right	\$125
ZARML	Flip-Up Arm Supports – 13"L pad (pr)	N/C		(18.5"L x 4.5"W – fits ZARML)	
ZALK	Arm Support Lock Kit (for Fixed Back) (pr)	\$86	ZPAML	Pivoting Arm Mount – Left (for ZFAPLL / ZTRL)	\$185
ZFAPLL	Long Molded Flat Armpad – Left (14"L x 3.25"W – fits ZARML)	\$90	ZPAMR	Pivoting Arm Mount – Right (for ZFAPLR / ZTRR)	\$185
ZFAPLR	Long Molded Flat Armpad – Right	\$90	ZPHGL	Pivoting Hand Grips (pr)	\$440
	(14"L x 3.25"W – fits ZARML)		ZAPSB*	Anterior Posture Support Bar	\$215
Z7FSL	Short Molded Flat Armpad – Left (10"L x 3"W – fits ZARMS)	\$90		* Enter back frame width to complete code (14" / 16" / 18") ¹	
Z7FSR	Short Molded Flat Armpad – Right (10″L x 3″W – fits ZARMS)	\$90	ZAPSB*	Anterior Posture Support Bar * Enter back frame width to complete code (20" / 22")1	\$250
			¹ Not compatible wi	th ZFL / ZFS / ZTR	

ZTAPLL	Tubular pad long 13"2 - Left (in lieu of molded pads, fits ZASWL)	N/A
ZTAPLR	Tubular pad long 13" ² - Right (in lieu of molded pads, fits ZASWL)	N/A
ZTAPSL	Tubular pad short 10° - Left (in lieu of molded pads, fits ZASWS)	N/A
ZTAPSR	Tubular pad short 10"2 - Right (in lieu of molded pads, fits ZASWS)	N/A

8. Foot and Leg Supports – Only ZFSL/R, ZFET2L/R, ZMFX8L/R, ZMFX12L/R are compatible with Zum Travel Kits (ZZUMAPTK)

ZFS		Foot Supports 16"- 20.5" (pr)	Std	ZRLSL	Universal Residual Limb Support – Left	\$215
ZFET	T2	Footrest Extension Tubes – 2" longer 1 (pr)	\$59	ZRLSR	Universal Residual Limb Support –	\$215
ZMF	-X8	MFX Foot Supports 8"-12" (pr)	\$115		Right	
ZMF	X12	MFX Foot Supports 12"-16.5" (pr)	\$115			
ZAD	M	Adj Angle/Depth Molded Footplates ² (pr)	\$215			
ZFET	T2A	Footrest Extension Tubes for Z7ADMFP (pr) – 2" longer	\$59	¹ Not compatible ² Not available fo	e with ZMFX8 / ZADM or ZMFX8	

9. Foot and Leg Support Accessories – ZFSLK, ZAHL, ZNCA, ZNCVS are NOT compatibe with Zum Travel Kits (ZZUMAPTK)

ZCS18	Calf Strap, Snap-On – Bodypoint Aeromesh – for 18"W base frame	\$150	ZAHL	Adjustable Heel Loops (pr) – not available for ZVS / ZFBT	\$69
ZCSIC18	Calf Strap, Snap-On – Infection Control – for 18"W base frame	\$140	ZNC	Neoprene Covers for Standard Footplates (Each Per Side)	\$51
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$100	ZNCABS16	Neoprene Covers for Adj Angle/Depth Molded Footplates (Each Per Side)	\$51



10. Positioning Pads -Nothing in Sec. 10 is compatible with Zum Travel Kits (ZZUMAPTK)

ZLTS1L	Lateral Thoracic Support w/ Small Pad ¹ – Left	\$365	ZLTS2L	Lateral Support w/ Large Depth- Adjustable Pad ¹ – Left	\$365
ZLTS1R	Lateral Thoracic Support w/ Small Pad ¹ – Right	\$365	ZLTS2R	Lateral Support w/ Large Depth- Adjustable Pad ¹ – Right	\$365
			ZLAEP	Lateral / Adductor 1.5" Extension Plate Quantity	\$27

¹ Depth-Adjustable Pad Sizes: Small – 5.5"x 3.5" / Large – 6.75" x 4.3". Pads can be adjusted medially up to 3" per side. For additional medial offset, specify ZLAEP extension plate.

9. Pelvic and Chest Belts - A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle	\$59 e)	ZCBM	Chest Belt ² – Medium (back widths 16"- 18", 1-piece w/ hook and loop closure)	\$42
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle	\$54)	ZCBL	Chest Belt ² – Large (back widths 18"- 22", 1-piece w/ hook and loop closure)	\$54
ZPBICM	Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	1 \$105	ZCBICM	Infection Control Chest Belt ^{1,2} – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	\$86
ZPBICL	Infection Control Pelvic Belt ¹ – Large (frame widths 18"-22", 2-piece w/ side-release buckle)		ZCBICL	·	\$96
ZPBBP_*	Pelvic Belt – Bodypoint Aeromesh * enter size to complete code (M / L)	\$285		W/ Side release buckley	
Recommended sizes for Pelvic belt ZCBBPM 16"-18 Bodypoint belts			8"W back fram	ne ZCBBPL 18"-22"W back frame	

¹Non-absorbent belts – recommended for infection-control environments

10. Other - ZPSK is not compatible with Zum Travel Kits (ZZUMAPTK)

ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$220	ZIHW	Whizard – winged urine deflector ²	\$74
ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$235	ZZUMAPCTLB		\$110
ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$260		Lanyard (set of 4)	
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$275			

¹Commode Pan Spacer Kits provide extra clearance for perianal care with pan in place. Reduces toilet bowl clearance height.

11. Credits and Omissions

Z905	Omit Commode Pan and	-\$20	Z907	Omit Standard Foot Supports	-\$45
	Mounting Brackets				



²Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

² Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.

SOME ASSEMBLY OF RAZ PRODUCTS IS REQUIRED.

12. Notes:

Specifications	Züm-AP		
Weight Capacity	350 lbs (160 kg)		
Overall Width – Frame Open	23.25" (59 cm)		
– Frame Folded	11" (28 cm)		
– Frame Folded	7.5" (19 cm) – (back assembly removed)		
Overall Length			
- with standard footrests	34" (86 cm)		
Seat Width	18" (46 cm)		
Clearance Height Under Seat			
– with pan brackets	17.25" (43 cm)		
– without pan brackets	17.75" (45 cm)		
Armrest Height Above Seat	8" (20 cm)		
Seat-to-Floor Height	21.5" (55 cm)		
Product Weight	44 lbs (20 kg)		
Frame Materials			
– base frame	Aluminum		
 back frame and footrests 	Stainless Steel		
– hardware	Stainless Steel		

