

Raz-AT ORDER FORM

EFFECTIVE FEBRUARY 2025 (CDN MSRP)

Date:__

| Dealer / Location: | | T: | _ F: |
|--------------------|----------|--------|------|
| PO# / Client: | Contact: | Email: | |

Raz-AT – Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height-adjustable stainless-steel frame, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. User weight capacity 400 lbs. See our <u>Evaluation and Medical Justification Form</u>.

| 1. Chair Codes | by size – <mark>Must Choose One</mark> | | | | | | | |
|------------------------------|--|---------|--|--|--|--|--|--|
| Z300S18 | Raz-AT w/ 18"W seat and base frame – Standard size, fits over most toilets ¹ | \$3,675 | | | | | | |
| Z300S* | Z300S* Non-Standard Width Raz-AT – Must add upcharge for custom width: see below * add custom width to complete code (16" / 20" / 22" / 24") – Adjusta-Back must be selected for all custom widths [Sec.6] | | | | | | | |
| | Z3SW16 Upcharge for custom 16"W seat and base frame ² – will not fit over most toilets | \$665 | | | | | | |
| | Z3SW20 Upcharge for custom 20"W seat and base frame ³ | \$665 | | | | | | |
| | Z3SW22 Upcharge for custom 22"W seat and base frame ⁴ | \$775 | | | | | | |
| | Z3SW24 Upcharge for custom 24"W seat and base frame ⁴ | \$875 | | | | | | |
| ³ 18"W seats can | its erface seats can be modified to fit. <mark>be selected but must also select Seat Support Adapter Kit ZSSKT [Sec.2]</mark> sco Foam Interface HD seat can be modified to fit. Seat Support Kit (ZSSKT) required for 24"W seat frame. | | | | | | | |
| Custom Seat / Ba | ase Frame Length and Seat Height | | | | | | | |
| Z3LNGSF3 | Upcharge for custom 3" Longer Seat Frame: required for seats \geq 21"D | \$665 | | | | | | |
| Z3LNGBF3 | Ipcharge for custom 3" Longer Base Frame: required for seats ≥ 22"D (typically used in \$665 onjunction with Z3SFL3) | | | | | | | |
| Z3LOWF5 | charge for custom Super-Low Frame (5" lower) – height adjustment range (molded seat): 16"-20" ¹ \$665 | | | | | | | |
| ¹ To lower height | adjustment range by 1" (15"-19"), select 4" caster upgrade option ZCUT4 (will not fit over most toilets) | | | | | | | |

2. Seats – Standard

| 2. Seats – Stanu | | | | | |
|-----------------------------------|--|----------|-------------------|--------------------------------------|-------|
| Molded Seats (thickness ~2.5"-3") | | | Visco-Foam Interf | ace Seats (thickness ~3.5 ″) | |
| ZMS17FRO | 18"W x 17"D – Front Access | Standard | ZVFI17BRI | 18"W x 17"D – Front Bridged | \$475 |
| ZMS17BRI | 18"W x 17"D – Front Bridged | N/C | ZVFI17FRO | 18"W x 17"D – Front Access | \$475 |
| ZMS20FRO | 18"W x 20"D – Front Access | \$150 | ZVFI17LTO | 18"W x 17"D – Left Access | \$475 |
| ZMS20BRI | 18"W x 20"D – Front Bridged | \$150 | ZVFI17RTO | 18"W x 17"D – Right Access | \$475 |
| Specialty Seats | | | ZVFI17RRO | 18"W x 17"D – Rear Access | \$475 |
| ZSOM | ZSOM Shower-Only Seat – Molded | \$440 | ZVFI17FAR | 18"W x 17"D – Front & Rear Access | \$475 |
| | 18″W x 16″-21″D | | ZVFI20BRI | 18″W x 20″D – Front Bridged | \$545 |
| ZEJRM | E&J Replica Seat – Molded – 18"W x 17"D | \$215 | ZVFI20FRO | 18"W x 20"D – Front Access | \$545 |
| | | ćopo | ZVFI20LTO | 18"W x 20"D – Left Access | \$545 |
| ZVFI23FRON | 23"W x 18"D HD – Visco Foam – Front Access Notched at rear to | \$930 | ZVFI20RTO | 18"W x 20"D – Right Access | \$545 |
| | maintain fore/aft adjustment | | ZVFI20RRO | 18"W x 20"D – Rear Access | \$545 |
| | range (IPAS) – 22"W frame | | ZVFI20FAR | 18"W x 20"D – Front & Rear Access | \$545 |
| h | | | | | |



3. Seats – Custom

| For custom seats, please spe Select Bridged if an access o | ecify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or F opening is not required. | Rear. | | | | | | |
|--|---|--------|--|--|--|--|--|--|
| ZCSEAT | Custom Seat (enter specs below) Select if both width and depth are $\leq 21''$ | \$1145 | | | | | | |
| ZCEJR | Custom E&J Replica Seat (enter specs below) Select if both width and depth are $\leq 21''$ | \$1305 | | | | | | |
| ZCSEATHD | Heavy Duty Custom Seats (enter specs below) Select if width and/or depth is $\geq 22''$ | \$1485 | | | | | | |
| Custom Specifications (for | above custom seats) | | | | | | | |
| Seat Frame Width: Custom Seat Width: Custom Seat Depth: Access Opening: Access Opening Width: Aperture Width: | Seat Frame Width: $ 16''/18''/20''/22''/24''$ Custom Seat Width: $ 16''/18''/20''/22''/24''$ (cannot be more than 2" wider than seat frame width)Custom Seat Depth: $ 15''/16''/17''/18''/19''/20''/21''/22''/23''/24'' [\geq 21'' requires Longer Frame (Z3SFL3)]Access Opening: FR - Front/LT - Left/RT - Right/RR - Rear/BR - BridgedAccess Opening Width: 2.5''/3.0''/3.5''/4.0''/4.5''/5.0'' (3.5''W is standard)$ | | | | | | | |
| Unless otherwise requested, the | width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats. | | | | | | | |
| Custom Seat width should not b | e narrower than frame width. | | | | | | | |
| Custom Seat can be up to 2" wid | er than seat frame width to provide additional thigh support without ordering wider chair frame. | | | | | | | |
| The Ischial Pelvic Alignment Syst | em (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats. | | | | | | | |
| Special modifications on Custom | Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc. | | | | | | | |

4. Seat Options

| ZSSKT | Seat Support Kit for Raz-AT (to support seats that are narrower | \$160 | ZGSB * Custom Seat Bracket Spacers – for \$325 adding dump angle (pr) |
|-------|--|-------|--|
| | than seat frames) | | Specify additional height at front of |
| | | | seat to complete code (0.5"/1.0"/1.5") |

5. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded, vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (85°/90°/95°) with cane-style posts. Back height is measured from top of molded seat to top of back upholstery. Hip Adductors [ZHP1/2] are not compatible with Fixed Back.

ZFB1818 Fixed Back – 18"W x 18"H^{1,2}

Standard

2 of 7



6. Adjusta-Back Frames

Includes height- and angle-adjustable arms with Armrest Lock Kit (ZALK), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (85°/90°/95°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. 24"H Adjusta-Back includes straight posts, on which the upholstery cannot be raised.

| Standard Width Adjust-Backs (18"W) | | | Non-Standard Width Adjusta-Backs | | | | | |
|------------------------------------|---|---------|----------------------------------|---|---|---------|--|--|
| ZAB1815 | Adjusta-Back – 15"H | \$275 | | * Complete Adjusta-Back order code by adding desired back frame width (16"/20"/22"/24") Note: 18" is the standard width. | | | | |
| ZAB1817 | Adjusta-Back – 17"H | \$275 | ZAB | | Adjusta-Back – 15"H non-std width | \$275 | | |
| ZAB1820 | Adjusta-Back – 20"H | \$355 | ZAB | - | Adjusta-Back – 17"H non-std width | \$275 | | |
| ZAB1824 | Adjusta-Back – 24"H | \$405 | ZAB | | Adjusta-Back – 20"H non-std width | \$355 | | |
| ZDTL | Dual Triggers – both on left side ¹ | N/C | ZAB | _24 * | Adjusta-Back – 24"H non-std width | \$405 | | |
| ZDTR | Dual Triggers – both on right side ¹ | N/C | ZRBT_ | _ 20* | Reclining Adjusta-Back ² 30° recline, | \$1,270 | | |
| ZRBT1820 | Reclining Adjusta-Back ² 30° recline, 20" back height with 2-Point | \$1,270 | | | 20" back height with 2-Point Reclining Arms (tilt is limited to 25°) | | | |
| | Reclining Arms (tilt is limited to 25°) | | ZRBT_ | _24* | Reclining Adjusta-Back ² 30° recline, | \$1,330 | | |
| ZRBT1824 | Reclining Adjusta-Back ² 30° recline, 24" back height with 2-Point Reclining Arms (tilt is limited to 25°) | \$1,330 | | | 24" back height with 2-Point Reclining Arms (tilt is limited to 25°) | | | |

¹ Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation. Available for all back widths. ²On reclining backs, dual triggers for recline are mounted on the left side and for tilt on the right side.

7. Custom Back Frame Widths Different From Seat Frames – Must choose back frame [Sec.5] (configurations below are not available for recliner backs)

| (j | clow are not available for recliner | | | | | | | | |
|---|--|-------|---|------------------|------------------|------------------|------------------|------------------|------|
| ZNB* | Custom Back – 2" Narrower than | \$380 | Seat Frame Width | 16″ ¹ | 18″ ² | 20″ ³ | 22″ ⁴ | 24″ ⁵ | |
| | Seat Frame | | Back Frame Width | 14″ | 16″ | 18″ | 20″ | 22″ | |
| | frame widths to complete code. e.g. ZNB2018 frame with an 18"W back frame | 3 | | | | | | | |
| ZNB* | Custom Back – 4" Narrower than | \$490 | Seat Frame Width | N/A | 18″ ¹ | 20″ ³ | 22″ ⁴ | 24″ ⁵ | |
| | Seat Frame – requires custom seat | | Back Frame Width | N/A | 14″ | 16″ | 18″ | 20″ | |
| | frame widths to complete code. s a 20″W seat frame with a 16″W back frame | | | | | | | | |
| ZWB* | Custom Back – 2″ Wider than Seat Frame | \$440 | Seat Frame Width | 16″ | 18″ | 20″ ³ | 22″ ⁶ | N/A | |
| | frame widths to complete code. | | Back Frame Width | 18″ | 20″ | 22″ | 24″ | N/A | |
| eg. ZWB1820 specifies an 18"W seat frame with a 20"W back frame | | | Consider also using Arm Spacer Kit (ZASK05/ZASK10) as another solution for providing more width between the arms. | | | | | | |
| 1 Deguires Custom See | 4 | | 1 | | | | | | |

¹ Requires Custom Seat

² Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

³ Custom Seat recommended. Also compatible with all 18"W seats, requires Seat Support Kit ZSSKT [Sec.2]

⁴Custom Seat recommended. Also compatible with modified ZVFI23FRO (18" seat depth with limited adjustability).

⁵Custom Seat recommended. Compatible with modified ZVFI23FRON (18" seat depth with limited adjustability, requires Seat Support Kit

(ZSSKT) [Sec.2] and modified ZVFI26FRO (18" seat depth with limited adjustability).

⁶ Custom Seat recommended. Also compatible with ZVFI23FRO (18" seat depth).



8. Back Support Options – These back supports replace the standard, tension-adjustable upholstery and are recommended for infection-control environments.

| ZSB14 | 14"W Symphony Back ¹ – recessed solid back | \$325 | | Harmony Back ² (contoured solid back | \$325 |
|-------|---|-------|-----------------------------|---|-------|
| ZSB16 | 16"W Symphony Back ¹ – recessed solid back | \$325 | | – 18"W x 13"H) Fits 17" or 20" back heights, 18"W back frames only | |
| ZSB18 | 18"W Symphony Back ¹ – recessed solid back | \$325 | | To W back frames only | |
| ZSB20 | 20"W Symphony Back ¹ – recessed solid back | \$440 | | | |
| ZSB22 | 22"W Symphony Back ¹ – recessed solid back | \$510 | | | |
| ZSB24 | 24"W Symphony Back ¹ – recessed solid back | \$535 | ¹ Not compatible | e with 15"H back – all Symphony Back panels are 13"H | |
| ZSB26 | 26"W Symphony Back ¹ – recessed solid back | \$555 | ² Not recommen | | |

9. Arm Supports and Accessories

| ZARML | Flip-Up Arm Supports – 13″L pad (pr) | Std | ZPAML | Pivoting Arm Mount – Left | \$185 |
|--|---|-------|--|---|-------|
| ZALK | Arm Support Lock Kit (for ZFB1818) (pr) | \$86 | | (For ZFAPLL/ZTRL) | |
| Z2PLT | 2-Point Flip-up Arm Supports ¹ (pr) (standard with reclining backs) | \$505 | ZPAMR | Pivoting Arm Mount – Right (for ZFAPLR/ ZTRR) | \$185 |
| ZFAPLL | Long Molded Flat Armpad ² – Left (14"L x 3.25"W) | \$90 | ZASK05 | Arm Support Spacer Kit – 1/2" per side (pr)(For Adjusta-backs) | \$69 |
| ZFAPLR | Long Molded Flat Armpad ² – Right | \$90 | ZASK10 | Arm Support Spacer Kit – 1″ per side (pr)(For Adjusta-backs) | \$130 |
| 7701 | (14″L x 3.25″W) | 4405 | ZPHGL | Pivoting Hand Grips (pr) | \$440 |
| ZTRL | Molded Arm Trough ² – Left (18.5″L x 4.5″W) | \$125 | ZAPSB* | Anterior Posture Support Bar * Enter back frame width to complete | \$215 |
| ZTRR | Molded Arm Trough ² – Right (18.5″L x 4.5″W) | \$125 | | code (14"/16"/18") ^{3,4} | |
| | | | ZAPSB* | Anterior Posture Support Bar * Enter back frame width to complete code (20"/22"/24") ^{3,4} | \$250 |
| ¹ For chairs with same back width and frame width only ² May require ZPAM for arm support to flip up and/or for lateral supports to swing away | | | ³ Not available with Z ⁴ Special length requi | FL / ZTR ired with ZASK05 (+1") and ZASK10 (+2"). | |



10. Foot and Leg Supports – Standard Foot Supports (ZFS) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support.

| ZFS | Foot Supports 16"-20.5" (pr) | Std | ZELRL | Elevating Leg Support ³ – Left | \$545 |
|----------|---|-------|-------|---|-------|
| ZFS16 | Foot Supports for 16"W frames | N/C | | (includes ZADML) | |
| 76670 | 16"-20.5" (pr) | ć c o | ZELRR | Elevating Leg Support ³ – Right | \$545 |
| ZFET2 | Footrest Extension Tubes – 2" longer ¹ (pr) | \$59 | | (includes ZADMR) | |
| ZMFX8 | MFX8 Foot Supports 8"-12" (pr) | \$115 | ZFOR | Forward Offset Foot Support Receiver (ea) Moves Foot support forward 2.5" | \$135 |
| ZMFX12 | MFX12 Foot Supports 12"-16.5" (pr) | \$115 | | to accommodate deeper seats (19" or deeper) Specify 1 or 2 | |
| ZADM | Angle/Depth Adj Molded Footplates ² (pr) | \$215 | ZRLSL | Universal Residual Limb Support – Left | \$215 |
| ZFET2A | Footrest Extension Tubes for ZADM (pr) (2" longer) | \$59 | ZRLSR | Universal Residual Limb Support – Right | \$215 |
| ZADABS16 | Angle/Depth Adj ABS Footplates (pr) For 16"W frames | \$215 | | ngn | |
| ZVS | V-Style Foot Support (1-piece) – for 16"/18"W base frames | \$325 | | | |
| ZVS* | V-Style Foot Support (1-piece) * Enter frame width to complete code (20" / 22" / 24") | \$360 | | | |
| ZFBT | Flip-Back Foot Support (1-piece) with calf pads | \$505 | | | |
| | vith ZMFX12 / ZADM / ZVS / ZFBAS vith Z300S16 or ZMFX8 | | | atible with frame widths ≥ 18″ and adjusted se ≤ 21″ w/ 3″ longer seat frames (Z3LNGSF3) | eat |

11. Foot and Leg Support Accessories

| | 5 | | | | |
|-------------------------------|--|-------------|----------|--|------|
| ZFSLK | Foot Support Lock (pr) | \$105 | ZAHL | Adjustable Heel Loops (pr) – not available for ZVS / ZFBT | \$69 |
| ZCS* | Snap-On Calf Strap – Bodypoint Aeromesh * Enter seat frame width to complete cod (16"/18"/20"/22"/24") | | ZNC | Neoprene Covers for Standard Footplates (Each Per Side) | \$51 |
| ZCSIC* | Calf Strap, Snap-On – Infection Control * Enter seat frame width to complete | \$140 ZNC16 | ZNC16 | Neoprene Covers for narrow molded footplates (Each Per Side) | \$51 |
| | code (16" / 18") | | ZNCA | Neoprene Covers for Angle/Depth Adj | \$51 |
| ZCSIC22 | Snap-On Infection Control Calf Strap – 20"/22"W frames | \$160 | ZNCABS16 | Molded Footplates (Each Per Side) Neoprene Covers for ABS Adj Angle/ | \$51 |
| ZCPL | Calf Panel – Bodypoint ¹ Fits frame widths 16" / 18" / 20" | \$200 | | Depth Footplates Z7ADABS (Each Per Side) | |
| ZHS1820 | H-Strap – for 18"/20"W frames | \$100 | ZLMOR | Lateral / Medial Offset Foot Support | \$90 |
| 21151020 | – hook and loop attachment | \$100 | | Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 | |
| ZHS2224 | H-Strap – for 22" / 24"W frames – hook and loop attachment | \$125 | | | |
| ¹ Not compatible w | vith ZMFX12 / ZFBT / ZVS | | | | |



12. Positioning Pads

| | J | | | | | |
|--------|--|-------|---|--|-------|--|
| ZHSP1 | Head Support Pad – Small (10" x 4") | Std | ZLGT1L | Leg Adductor w/ Small Pad ¹ – Left | \$365 | |
| ZHSP2 | Head Support Pad – Large (12" x 6") | \$59 | ZLGT1R | Leg Adductor w/ Small Pad ¹ – Right | \$365 | |
| ZOHIP | Lateral Offset Head Support Interface Plate – up to 22° of adjustment | \$86 | ZLGT2L | Leg Adductor w/ Large Pad ¹ – Left | \$365 | |
| ZLTS1L | | \$365 | ZLGT2R | Leg Adductor w/ Large Pad ¹ – Right | \$365 | |
| ZLIJIL | Lateral Thoracic Support w/ Small Pad ^{1,2} – Left | 200 | ZHP1L | Hip Adductor w/ Small Pad ^{1,3} – Left | \$365 | |
| ZLTS1R | R Lateral Thoracic Support w/ Small Pad ^{1,2} | \$365 | ZHP1R | Hip Adductor w/ Small Pad ^{1,3} – Right | \$365 | |
| | – Right | | ZHP2L | Hip Adductor w/ Large Pad ^{1,3} – Left | \$365 | |
| ZLTS2L | Lateral Thoracic Support w/ Large Pad ^{1,2} | \$365 | ZHP2R | Hip Adductor w/ Large Pad ^{1,3} – Right | \$365 | |
| | – Left | | ZMTST1 | Medial Thigh Support – Small | \$495 | |
| ZLTS2R | Lateral Thoracic Support w/ Large Pad ^{1,2} – Right | \$365 | ZMTST2 | Medial Thigh Support – Large | \$495 | |
| ZLAEP | Lateral/Adductor Extension Plate | \$27 | ¹ Adjustable Pad Sizes: Small – 5.5" x 3.5" / Large – 6.75" x 4.3". ² Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side For additional medial offset, specify extension plate. ³ Fits Adjusta-Back only | | | |

13. Pelvic and Chest Belts – A chest belt should not be used without a pelvic belt

| Pelvic Belts | | | | | Chest Belts | | | | |
|--------------------------|--|--------|-----------------------|---------|--|---|------------|------------------------|-------|
| ZPBS | Pelvic Belt – Small (frame widths 14"-16", 2-piece w/ side-release buckle) | | \$48 | ZCBS | Chest Belt ² – Small (back widths 14"-16 1-piece w/ hook and loop closure) | | | \$37 | |
| ZPBM | Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle) | | \$59 | ZCBM | Chest Belt ² – Medium (back widths 16"- 18", 1-piece w/ hook and loop closure) | | | \$42 | |
| ZPBL | Pelvic Belt – Large (fra 18"-22", 2-piece w/ sid | | uckle) | \$54 | ZCBL | Chest Belt ² – Lar 22", 1-piece w/ h | J . | | \$54 |
| ZPBXL | Pelvic Belt – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle) | | \$74 | ZCBXL | Chest Belt – XLarge ² (back frame widths 22"-30", 1-piece w/ hook and loop | | ne widths | \$69 | |
| ZPBICS | Infection Control Pelvic Belt ¹ – Sm | | nall | \$105 | | closure) | | 500 | |
| | (frame widths 14"-16", 2-piece w/ side-release buckle) | | | Z | ZCBICS | Infection Control Chest Belt ^{1,2} – Small (back widths 14-16", 1-piece | | – Small | \$86 |
| ZPBICM | Infection Control Pelvic Belt ¹ – Medium (frame widths 16"-18", 2-piece w/ side-release buckle) | | edium | m \$105 | | w/ side-release b | ouckle) | | |
| | | | | | ZCBICM | Infection Contro | | | |
| ZPBICL | Infection Control Pelvic Belt ¹ – Large | | \$105 | | um (back widths 16"-18", 1-piece w/ side-release buckle) | | ce | | |
| | (frame widths 18"-22", 2-piece w/ side-release buckle) | | | ZCBICL | Infection Control Chest Belt ^{1,2} – Large (back widths 18"-22", 1-piece | | | \$96 | |
| ZPBICXL | Infection Control Pelvic Belt ¹ – XLarge | | | \$120 | | w/ side-release buckle) | | | |
| | (seat frame widths 22"-30", 2-piece w/ side-release buckle) | | | ZCBICXL | Infection Control Chest Belt ^{1,2} – XLarge (back frame widths 22"-30", 1-piece w/ side-release buckle) | | | \$120 | |
| ZPBBP_* | Pelvic Belt – Bodypoint Aeromesh \$2 * enter size to complete code (M / L) | | \$285 | | | | | | |
| ZPBBPXL | Pelvic Belt – Bodypoint | | | \$215 | ZCBBP_* | Chest Belt – Body * enter size to co | • | | \$295 |
| Recommended sizes for | d Pelvic belt | ZPBBPM | 16″/18"W seat fram | | ZPBBPL | 18"-22"W seat frame | ZPBBPXL | 22"-26"W seat frame | |
| Bodypoint belts | Chest belt | ZCBBPM | 16"/18"V back frar | - | ZCBBPL | 18"/20"W back frame | ZCBBPXL | 22″-26"W | |

²Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.



7 of 7

| 14. Other | | | , | KA4 | Z-AT 09-01-25-0 |
|---|--|-------|---|---|-----------------|
| ZDCLT | Directional Caster Lock – fits on one rear | \$185 | ZIHW | Whizard – winged urine deflector ² | \$74 |
| | caster (facilitates steering) | | ZBUDBH* | Bard Urine Drainage Bag Holder | \$180 |
| ZATTRR | Rear anti-Tippers for 5" casters (pr) | \$380 | | * Enter size to complete code | |
| ZCUT4 | 4" Caster Upgrade (4x) | \$135 | | (16"/18"/20"/22"/24") | |
| ZCUT6 | 6" Caster Upgrade (4x) | \$385 | ZSPKT22 | AT Self-Propel Kit with 22" wheels | \$1695 |
| ZPSK1 | Commode Pan Spacer Kit ¹ – 1″ lower | \$220 | – incl anti-tippers (tilt limited to 2) | |) |
| ZPSK2 | Commode Pan Spacer Kit ¹ – 2" lower | \$235 | ZSPKT24 | AT Self-Propel Kit with 24" wheels – incl anti-tippers (tilt limited to 25°) | \$1695 |
| ZPSK3 | Commode Pan Spacer Kit ¹ – 3" lower | \$260 | | | |
| ZPSK4 | Commode Pan Spacer Kit ¹ – 4" lower | \$275 | | | |
| ZTLS Tilt Limit Spacers – 1/2" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for addi- tional tilt limiting. Quantity | | N/C | | bacer Kits provide extra clearance for self- ces toilet bowl clearance height. | care with |
| | | | ² Has patent-pendi medial aspect of th | ng lateral wings that are positioned unde thighs to keep the Whizard in place. | r the |

15. Credits and Omissions

| Z905 | Omit Commode Pan and Mounting Brackets | -\$20 | Z911 | Omit Molded Head Support Pad | -\$20 |
|------|---|-------|------|-------------------------------------|-------|
| Z907 | Omit Foot Supports - if no other Foot / Leg Support Option is selected | -\$45 | Z912 | Omit Complete Head Support Assembly | -\$50 |

SOME ASSEMBLY OF RAZ PRODUCTS IS REQUIRED.

| Specifications | Raz-AT (Attendant Tilt) | | | | |
|---|---|--|--|--|--|
| Weight Capacity | 400 lbs | | | | |
| Overal Width | Seat width plus 4.5" with Standard Arm Pads | | | | |
| Overall Length (with foot supports and headrest adj | usted to middle positions) | | | | |
| - at 0° tilt and standard foot supports | 43″ | | | | |
| - at 40° | 54" | | | | |
| Seat Depth | 16" or 21" depending on seat/cushion length | | | | |
| Toilet Clearance Height (under seat) | | | | | |
| - with pan brackets | 17"-21" (1" increments) | | | | |
| - without pan brackets | 17.75"-21.75 (1" increments) | | | | |
| Toilet Clearance Width | 15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25" | | | | |
| Arm Support Height - above molded seat (1/2" incre | ments) | | | | |
| - 18"H Fixed Back - standard | 8.5″ | | | | |
| - 15″H Adjusta-Back | 5.5"-8.5" (1/2" increments) | | | | |
| - 17"H Adjusta-Back | 6.5"-9.5" (1/2" increments) | | | | |
| - 20"H Adjusta-Back | 7.5"-10.5" (1/2" increments) | | | | |
| - 24"H Adjusta-Back | 7.5"-10.5" (1/2" increments) | | | | |
| - 20"H Reclining Adjusta-Back | 7.5"-10.5" (1/2" increments) | | | | |
| - 24"H Reclining Adjusta-Back | 7.5"-10.5" (1/2" increments) | | | | |
| Seat-to-Floor Height (with molded seat) | 21"-25" in 1" increments | | | | |
| Product Weight (with 18"W seat) | 61 lbs | | | | |
| Frame and Hardware Material | Stainless Steel (lifetime warranty against corrosion) | | | | |

