

RAZE ORDER FORM RAISED TOILET SEAT / TRANSFER TUB BENCH


CANADIAN **EFFECTIVE APRIL 2025**

Date:		


Dealer / Location:		T:	F:
PO# / Client:	Contact:	Email:	

Raze modular frame for over-toilet use / tub transfer bench when combined with the optional transfer surface. User weight capacity 352 lbs.

Chair Code	Seat Frame Width	MSRP
Z4441	Raze Frame with tall leg set – for over-toilet / shower stall use	\$1,435
Z4442	Raze Frame with short leg set – for tub use MUST SELECT TRANSFER SURFACE BELOW	\$1,485
Transfer Surface Optional for sh	ces (for tub and shower stall use) – MUST select a Transfer Surface for use in tub (Z4442). ower stall use.	
ZRZTSN	Transfer Surface – Narrow (8.5" x 16"D) [Suitable for most tubs]	\$645
ZRZTSW	Transfer Surface – Wide (12.5" x 16"D)	\$765







Z4442 with optional ZRB1815 and ZRZTSN

2. Optional Components

ZRZBS1815	15"H Back Support including swing-away arm supports (Recommened for tub use)	\$445	
ZCHL	Removable Cane Handle – Left (Required when Chest Belt is ordered)	\$75	
ZCHR	Removable Cane Handle – Right (Required when Chest Belt is ordered)	\$75	
Additional Leg Sets (for use over toilet, in shower stall and in tub			
ZRZSL	Short Leg Set – for tub use (when Z4441 is ordered)	\$175	
ZRZTL	Tall Leg Set – for toilet use / shower stall use with Transfer Surface (when Z4442 is ordered)	\$135	



3. Seats – Standard

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interf	ace Seats (thickness ~3.5 ″)	
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18″W x 17″D – Front Bridged	\$475
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$475
	5		ZVFI17LTO	18"W x 17"D – Left Access	\$475
ZMS20FRO	18"W x 20"D – Front Access	\$150	ZVFI17RTO	18"W x 17"D – Right Access	\$475
ZMS20BRI	18"W x 20"D – Front Bridged	\$150	ZVFI17RRO	18"W x 17"D – Rear Access	\$475
Specialty Seats			ZVFI20BRI	18″W x 20″D – Front Bridged	\$545
ZSOM	Shower-Only Seat – Molded	\$440	ZVFI20FRO	18"W x 20"D – Front Access	\$545
250111	18"W x 16"-21"D	ψ Π O	ZVFI20LTO	18"W x 20"D – Left Access	\$545
7F JRM	E&J Replica Seat – Molded	\$215	ZVFI20RTO	18"W x 20"D – Right Access	\$545
225.001	– 18"W x 17"D	<i>7213</i>	ZVFI20RRO	18"W x 20"D – Rear Access	\$545

4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required.

ZCSEAT	Custom Seat (enter specs below) Select if both width and depth are $\leq 21''$	\$1145
ZCEJR	Custom E&J Replica Seat (enter specs below) Select if both width and depth are $\leq 21''$	\$1305
Custom Specifications (for a	above custom seats)	
Seat Frame Width: Custom Seat Width: Custom Seat Depth: Access Opening: Access Opening Width: Aperture Width:	 18" 18"/20" 17"/18"/19"/20" FR - Front / LT - Left / RT - Right / RR - Rear / BR - Bridged 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard) 3.0"-8.0" (6.5" recommended) - enter "NA" if aperture is not required 	
Unless otherwise requested, the	width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.	
Custom Seat width should not be	e narrower than frame width.	
Custom Seat can be up to 2" wide	er than seat frame width to provide additional thigh support without ordering wider chair frame.	
The Ischial Pelvic Alignment Syste	em (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.	
Spacial modifications on Custom	South include: log longth discremency coff spots prossure relief areas, contouring atc	

Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

5. Arm Supports - Wrap-around, swing-away, removable

ZASWL Arm Supports (p) Wraparound – Long	Standard
11 4		Standard
ZASWS Arm Supports (p) Wraparound – Short	N/C
Arm Support Pads		
ZARMLL Molded pad long	13" – Left (fits ZASWL)	Standard
ZARMLR Molded pad long	13" – Right (fits ZASWL)	Standard
ZARMSL Molded pad shore	t ¹ 10″ – Left (fits ZASWS)	N/C
ZARMSR Molded pad shore	t¹ 10″ – Right (fits ZASWS)	N/C
ZTAPLL Tubular pad long	13" – Left (in lieu of molded pads, fits ZASWL)	N/C
ZTAPLR Tubular pad long	13" – Right(in lieu of molded pads, fits ZASWL)	N/C
ZTAPSL Tubular pad shor	t 10" – Left (in lieu of molded pads, fits ZASWS)	N/C
ZTAPSR Tubular pad shor	t 10" – Right (in lieu of molded pads, fits ZASWS)	N/C



3 of 4

ZFAPLL	Molded Flat Armpad Long 14"L x 3.25"W – Left (fits ZASWL)	\$90
ZFAPLR	Molded Flat Armpad Long 14"L x 3.25"W – Right (fits ZASWL)	\$90
ZFAPSL	Molded Flat Armpad Short 10"L x 3"W – Left (fits ARMS)	\$90
ZFAPSR	Molded Flat Armpad Short 10"L x 3"W – Right (fits ARMS)	\$90
ZTRL	Molded Arm Trough 18.5"L x 4.5"W – Left (fits ZASWL)	\$125
ZTRR	Molded Arm Trough 18.5″L x 4.5″W – Right (fits ZASWL)	\$125
ZPAML	Pivoting Arm Mount ² – Left (for ZFAPLL / ZTRL)	\$185
ZPAMR	Pivoting Arm Mount ² – Right (for ZFAPLR / ZTRR)	\$185
ZAPSB18	Anterior Posture Support Bar (for ZARML only)	\$215
¹ Not compatik	ble with ZFAP / ZTR / ZPAM	

² Recommended when Lateral Thoracic Supports (LTL/S) and Molded Arm Troughs (ZTR) are both ordered

6. Pelvic and Chest Belts – A chest belt should not be used without	a pelvic belt
	a pertie sere

Pelvic Belts					Chest Belts		
ZPBM	Pelvic Belt – M 16″-18″, 2-pie			• • •	ZCBM	Chest Belt ² – Medium (back widths 16"- 18", 1-piece w/ hook and loop closure)	\$42
ZPBL		Pelvic Belt – Large (frame widths \$54 8"-22", 2-piece w/ side-release buckle)			ZCBL	Chest Belt ² – Large (back widths 18"- 22", 1-piece w/ hook and loop closure)	\$54
ZPBICM	Infection Cor (frame width: w/ side-relea:	s 16″-18″, 2-		ium \$105	ZCBICM	Infection Control Chest Belt ^{1,2} – Medi- um (back widths 16"-18", 1-piece w/ side-release buckle)	\$86
ZPBICL	Infection Cor (frame width: w/ side-relea:	s 18″-22″, 2-		e \$105	ZCBICL	Infection Control Chest Belt ^{1,2} – Large (back widths 18"-22", 1-piece w/ side-release buckle)	\$96
ZPBBP_*	Pelvic Belt – E * enter size to			\$285)	ZCBBP_*	Chest Belt – Bodypoint Aeromesh * enter size to complete code (M / L / XL)	\$295
Recommended sizes for	d Chest belt	ZPBBPM	16″/18"W	seat frame	ZPBBPL	18"-22"W seat frame	
Bodypoint belts	Chest Delt	ZCBBPM	16"/18"W	back frame		18"/20"W back frame	
¹ Non-absorbent	belts – recommer	nded for infecti	on-control en	vironments			

¹Non-absorbent belts – recommended for infection-control environments ²Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.

7. Positioning Pads

	ic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. medial offset, select ZLAEP extension plate.	
ZLTSWA1L	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Left	\$365
ZLTSWA1R	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Right	\$365
ZLTSWA2L	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Left	\$365
ZLTSWA2R	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Right	\$365



8. Urine Management

ZIHW	Whizard – winged urine deflector – Has patented lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$74
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$165
ZSC17	Splash Curtain Fixed Style for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$330
ZSC20	Splash Curtain Fixed Style for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$375
ZSCC	Custom Splash Curtain Fixed Style – Raises seat 1/2"	\$405
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$375
ZCPAN	Commode Pan (with brackets)	N/C
¹ For Left, Rig	ght, Front, and Rear Open, call for order codes and pricing.	

Transfer Tub Bench Raised Toilet Seat Transfer Tub Bench Specifications Raised Toilet Seat w/ Optional Back Support w/ Optional Back Support Weight Capacity 352 lbs **Overall Width** 20"-23.75" 23.25"-23.75" - w/ Narrow Transfer Surface 29.75"-31.25" (8.5") - w/ Wide Transfer Surface 33.75"-35.25" (12.5") Overall Lenght 20"-24.25" (front to rear) 15″ 15″ Back support height N/A N/A above seat Back support Width N/A 18″ N/A 18″ Arm support height 9" (w/ molded seat) above seat 19.5″ Seat frame weight Seat depth 17" or 20" (middle of IPAS range)

Specifications	Raised Toilet Seat	Raised Toilet Seat w/ Optional Back Support	Transfer Tub Bench	Transfer Tub Bench w/ Optional Back Support
Clearance height under seat- 0.5" increments	16.5"-19.5" (w/ molded seat) * without brackets +0.75"		N/A	
Under-seat clearance width	11.25″-15.5″		N/A	
Overall Lenght (front to rear)	20"-24.25"			
Seat-to-floor height	20.25"-23.25" 0.5" increments		N/A	
Transfer surface-to-floor height	N/A		21"-24" (to floor outside tub)	
Tub wall clearanceheight	N/A		16.75"-19.75" (to floor outside tub)	
Product weight	19 lbs	26 lbs	26.5 lbs	33.5 lbs
Frame and hardware material	Aluminum			
¹ Lifetime warranty against corrosion				

